

Vitalograph Launch Powerful New In2itive e-Diary for Global Clinical Trials

Next generation eCOA device with integrated respiratory assessment



Vitalograph are delighted to announce that their most powerful ever In2itive e-Diary is now available for global clinical trials, offering the ideal solution for respiratory endpoint studies.

The new In2itive™ e-Diary is a robust handheld medical device with touchscreen for eCOA data gathering. Using Vitalograph's proven Fleisch Pneumotachograph flowhead, it provides accurate measurement of respiratory endpoints in clinic or for home monitoring.

Features

- Easy to use with large, high resolution, responsive touch screen.
- Live spirometry feedback promotes high quality, reliable data.
- Integrated training modules provide technique reminders supporting compliance and data quality.
- Alerts, reminders, and workflows are fully customizable to fit study protocols.
- Pre validated for a range of e-PROs to cover respiratory endpoints and health status.
- Date/time stamps ensuring data integrity.
- Secure end-to-end encrypted data transmission.

Fully integrated wireless communications (cellular/mobile) allow automatic secure transmission of data to the study web portal. Uploaded data is available for instant review, giving remote insight into site and subject performance. The In2itive e-Diary is a fully validated class II (USA) / IIa (EU) medical device, meeting or exceeding all technical requirements for cyber-security and data protection.

The new device meets the performance requirements of ISO 23747 (peak flow), ISO 26782 (volumes) and adheres to guidance in the ATS/ERS joint 2019 technical statement on the standardization of spirometry.

About Vitalograph

Vitalograph is a global leader in respiratory diagnostic devices, clinical trial services and medical equipment servicing. With a pioneering heritage of excellence spanning half a century, Vitalograph continues to make valuable contributions to effective medical care and enhanced quality of life.

[Watch the In2itive e-Diary video here](#)

Press contact | Helen Venn | helen.venn@vitalograph.co.uk